

Bosentan Hydrate

Catalog No: #S6265


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bosentan Hydrate
Brief Description	Inhibitors
Purification	/
Target Name	Endothelin Receptor antagonist
Calculated MW	569.63
Formulation	C ₂₇ H ₂₉ N ₅ O ₆ S.H ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: O.COC1=CC=CC=C1OC1=C(OCCO)N=C(N=C1NS(=O)(=O)C1=CC=C(C=C1)C(C)(C)C)C1=NC=CC=N1
 [c:4,6,10,16,18,27,29,38,40,t:2,25,36]
 Pathways: GPCR/G Protein Receptor: ET-A, ET-B
 Boiling pt: Melting pt: Solubility: DMSO: 100 mg/mL (175.55 mM) Ethanol: 2 mg/mL (3.51 mM) Water: 0.001 mg/mL (0.0 mM)
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only