

## Trichostatin A (TSA)

Catalog No: #S6270


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Trichostatin A (TSA)
Brief Description	Inhibitors
Purification	99.02%
Target Name	HDAC inhibitor
Calculated MW	302.4
Formulation	C17H22N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer SMILES: C[C@H](/C=C(\C)/C=C/C(=O)NO)C(=O)C1=CC=C(C=C1)N(C)C Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair; NF-kb Receptor: HDAC

Boling pt: Melting pt: Solubility: DMSO: 23 mg/mL (76.05 mM) Ethanol:

## References

1. Vigushin DM, et al. Clin Y Res, 2001, 7(4), 971-976.

Note: This product is for in vitro research use only