

## Tipifarnib


Catalog No: #S6271

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tipifarnib
Brief Description	Inhibitors
Purification	99.00%
Target Name	Transferase inhibitor
Calculated MW	489.4
Formulation	C27H22Cl2N4O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer SMILES: CN1C=NC=C1[C@@](C2=CC=C(C=C2)Cl)(C3=CC4=C(C=C3)N(C(=O)C=C4C5=CC(=CC=C5)Cl)C)N Pathways: Metabolism Receptor: Farnesyl transferase (FTase) Boiling pt: Melting pt: Solubility: DMSO: 14 mg/mL (28.6 mM) Ethanol:

## References

1. Margolin KA, et al. Clin Y Res. 2012, 18(4), 1129-1137.

Note: This product is for in vitro research use only