

Enzastaurin (LY317615)

Catalog No: #S6280


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Enzastaurin (LY317615)
Brief Description	Inhibitors
Purification	98.00%
Target Name	PKC inhibitor
Calculated MW	515.61
Formulation	C32H29N5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CN1C=C(C2=CC=CC=C12)C1=C(C(=O)NC1=O)C1=CN(C2CCN(CC3=NC=CC=C3)CC2)C2=CC=CC=C12
 |c:2,6,30,32,40,t:4,8,12,20,28,38,42|Pathways: Cytoskeletal Signaling Receptor: PKCα; PKCβ; PKCγ; PKCε
 Boling pt: Melting pt: Solubility: DMSO: 30 mg/mL (58.18 mM) Ethanol:

References

1. Graff JR, et al. Cancer Res, 2005, 65(16), 7462-7469.

Note: This product is for in vitro research use only