

Tariquidar

Catalog No: #S6287


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tariquidar
Brief Description	Inhibitors
Purification	99.00%
Target Name	P-gp inhibitor
Calculated MW	646.73
Formulation	C38H38N4O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: COC1=CC2=C(CN(CCC3=CC=C(NC(=O)C4=C(NC(=O)C5=CC6=CC=CC=C6N=C5)C=C(OC)C(OC)=C4)C=C3)CC2)C=C1OC

|c:17,26,28,31,40,43,49;t:2,4,10,12,22,24,34|Pathways: Membrane transporter/Ion channel; Neuroscience Receptor: P-gp?(P-glycoprotein) Boling

pt: Melting pt: Solubility: DMSO: 52 mg/mL (80.4 mM) Ethanol:

References

1. Martin C, et al. Br J Pharmacol, 1999, 128(2), 403-411.

Note: This product is for in vitro research use only