

Dabigatran (BIBR 953)

Catalog No: #S6295


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dabigatran (BIBR 953)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Thrombin inhibitor
Calculated MW	471.51
Formulation	C25H25N7O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: CN1C2=C(C=C(C=C2)C(=O)N(CCC(=O)O)C3=CC=CC=N3)N=C1CNC4=CC=C(C=C4)C(=N)N Pathways: Others Receptor: Thrombin
 Boiling pt: Melting pt: Solubility: DMSO: 0.5mg/ml(1.06 mM); Water:

References

1. Huel NH, et al. J Med Chem. 2002, 45(9), 1757-1766.

Note: This product is for in vitro research use only