

Lonafarnib

Catalog No: #S6302


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lonafarnib
Brief Description	Inhibitors
Purification	98.00%
Target Name	Raf inhibitor; Ras inhibitor
Calculated MW	638.82
Formulation	C27H31Br2ClN4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: NC(=O)N1CCC(CC(=O)N2CCC(CC2)[C@H]2C3=C(CCC4=C2C(Br)=CC(Cl)=C4)C=C(Br)C=N3)CC1
 [r,c:18,22,26,29,35,t:32|Pathways:MAPKReceptor:K-ras-4B; H-ras; N-ras
 Boling pt: Melting pt: Solubility: DMSO: 127 mg/mL (198.8 mM) Ethanol: 127 mg/mL (198.8 mM) Water:

References

1. Liu M, et al. Cancer Res, 1998, 58(21), 4947-4956.

Note: This product is for in vitro research use only