

Lubiprostone

Catalog No: #S6307


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lubiprostone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	390.46
Formulation	C ₂₀ H ₃₂ F ₂ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: CCCCC([C@]1(CC[C@H]2[C@H](O1)CC(=O)[C@@H]2CCCCC(=O)O)O)(F)F

Pathways: Others Receptor: CIC-2 chloride channels

Boiling pt: Melting pt: Solubility: DMSO: 78 mg/mL (199.76 mM) Ethanol: 78 mg/mL (199.76 mM) Water

References

1. Bijvelds MJ, et al. Gastroenterology, 2009, 137(3), 976-985.

Note: This product is for in vitro research use only