

Stattic

Catalog No: #S6308


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Stattic
Brief Description	Inhibitors
Purification	99.00%
Target Name	STAT inhibitor
Calculated MW	211.19
Formulation	C ₈ H ₅ NO ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: [O-][N+](=O)C1=CC=C(C=C1)C=C2C(=O)N(S(=O)(=O)C2)C=C1 Pathways: JAK/STAT signaling; Stem Cells Receptor: STAT3 Boiling pt: Melting pt: Solubility: DMSO: 42 mg/mL (198.87 mM) Appearance: Orange Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Schust J, et al. Chem Biol, 2006, 13(11), 1235-1242.

Note: This product is for in vitro research use only