


Description

Product Name	GW3965 HCl
Brief Description	Inhibitors
Purification	98.00%
Target Name	Liver X Receptor agonist
Calculated MW	618.51
Formulation	C33H31ClF3NO3.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:Cl.O=C(O)CC1=CC=CC(OCCCN(CC(C2=CC=CC=C2)C2=CC=CC=C2)CC2=CC=CC(=C2Cl)C(F)(F)F)=C1

[c:6,18,20,25,27,33,35,42,t:4,16,23,31|Pathways:OthersReceptor:hLXRα; hLXRβ; LXRα/SRC1?LiSABoling pt:Melting pt:Solubility:DMSO:16 mg/mL warmed (25.86 mM)Ethanol:

References

1. Collins JL, et al. J Med Chem. 2002, 45(10), 1963-1966.

Note: This product is for in vitro research use only