

Tie2 kinase inhibitor


Catalog No: #S6335

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tie2 kinase inhibitor
Brief Description	Inhibitors
Purification	98.00%
Target Name	Tie-2 inhibitor
Calculated MW	439.53
Formulation	C ₂₆ H ₂₁ N ₃ O ₂ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: COC1=CC2=C(C=C1)C=C(C=C2)C1=C(N=C(N1)C1=CC=C(C=C1)S(C)=O)C1=CC=NC=C1

|c:4,6,9,11,16,22,24,32,34,t:2,14,20,30|Pathways: Tyrosine Kinase/Adaptors Receptor: Tie-2 Boiling pt: Melting pt: Solubility: DMSO: 35 mg/mL (79.63 mM) Ethanol:

Note: This product is for in vitro research use only