


## Description

Product Name	KX2-391
Brief Description	Inhibitors
Purification	99.43%
Target Name	Src inhibitor
Calculated MW	431.53
Formulation	C26H29N3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer SMILES: O=C(CC1=CC=C(C=N1)C1=CC=C(OCCN2CCOCC2)C=C1)NCC1=CC=CC=C1

|c:5,7,24,31,33,t:3,10,12,29|Pathways: Angiogenesis Receptor: Src (Hep 3B); Src (Hep G2); Src (HuH7); Src (PLC/PRF/5) Boling pt: Melting

pt: Solubility: DMSO: 86 mg/mL (199.29 mM) Ethanol:

## References

1. Lau GM, et al, Dig Dis Sci, 2009, 54(7), 1465-1474.

Note: This product is for in vitro research use only