

# (R)-Nepicastat HCl


Catalog No: #S6356

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	(R)-Nepicastat HCl
Brief Description	Inhibitors
Purification	/
Target Name	Hydroxylase inhibitor
Calculated MW	331.81
Formulation	C <sub>14</sub> H <sub>15</sub> F <sub>2</sub> N <sub>3</sub> S.HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system SMILES: Cl.NCC1=CNC(=S)N1[C@@H]1CCC2=C(F)C=C(F)C=C2C1

[r,c:12,18,t:2,15] Pathways: Metabolism Receptor: Dopamine β-hydroxylase; Dopamine β-hydroxylase Boiling pt: Melting pt: Solubility: DMSO: 66 mg/mL (198.9 mM) Ethanol:

Note: This product is for in vitro research use only