

## 1-Hexadecanol


Catalog No: #S6359

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	1-Hexadecanol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	242.44
Formulation	C16H34O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:CCCCCCCCCCCCCCCCOPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:68 mg/mL (280.48 mM);Ethanol:68 mg/mL (280.48 mM);Water:68 mg/mL (280.48 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only