

3-Indolebutyric acid (IBA)


Catalog No: #S6361

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3-Indolebutyric acid (IBA)
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	203.24
Formulation	C ₁₂ H ₁₃ NO ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: OC(=O)CCCC1=CNC2=C1C=CC=C2 | c:9,12,14,t:6 | Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 41 mg/mL (201.73 mM) Ethanol: 41 mg/mL (201.73 mM) Water:

Note: This product is for in vitro research use only