

# 5-hydroxymethyl Tolterodine (PNU 200577, 5-HMT, 5-HM)


Catalog No: #S6364

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	5-hydroxymethyl Tolterodine (PNU 200577, 5-HMT, 5-HM)
Brief Description	Inhibitors
Purification	/
Target Name	AChR antagonist
Calculated MW	341.49
Formulation	C22H31NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system SMILES: CC(C)N(CC[C@H](C1=CC=CC=C1)C1=C(O)C=CC(CO)=C1)C(C)C

[r,c:9,11,14,17,21,t:7] Pathways: Neuroscience Receptor: mAChR Boling pt: Melting pt: Solubility: DMSO: 68 mg/mL (199.12 mM) Ethanol: 68 mg/mL (199.12 mM) Water:

Note: This product is for in vitro research use only