

Abscisic Acid (Dormin)


Catalog No: #S6368

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Abscisic Acid (Dormin)
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	264.32
Formulation	C15H20O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC(=C)C(C)C(=O)C1(C)C=C(O)O | r,c:7|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:52 mg/mL (196.73 mM)Ethanol:52 mg/mL (196.73 mM)Water:

Note: This product is for in vitro research use only