

## Alibendol

Catalog No: #S6375


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Alibendol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	251.28
Formulation	C13H17NO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:COC1=C(O)C(=CC(CC=C)=C1)C(=O)NCCOPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 50 mg/mL (198.98 mM); DMSO: 50 mg/mL (198.98 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Sheng, L., H. Chen, and Y. Li,et al. J Chromatogr B Analyt Technol Biomed Life Sci, 2007. 854(1-2): p. 99-103.

Note: This product is for in vitro research use only