

Allopurinol Sodium


Catalog No: #S6376

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Allopurinol Sodium
Brief Description	Inhibitors
Purification	/
Target Name	ROS inhibitor
Calculated MW	159.1
Formulation	C ₅ H ₄ N ₄ NaO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: [Na+].O=C1NC=NC2=C1C=NN2 [c:3,5,8]
Pathways: Immunology/Inflammation
Receptor: XOB
Boiling pt: Melting pt: Solubility: DMSO: 32 mg/mL (201.13 mM) Ethanol: 2 mg/mL (12.57 mM) Water: 32 mg/mL (201.13 mM)
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only