

## Aloperine

Catalog No: #S6377


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Aloperine
Brief Description	Inhibitors
Purification	99.66%
Target Name	Others
Calculated MW	232.36
Formulation	C15H24N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: C1CCN2C[C@H]3C[C@@H]([C@H]2C1)C=C4[C@H]3NCCC4  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: Ethanol: 46 mg/mL (197.96 mM); DMSO: 2 mg/mL (8.6 mM)  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Yuan XY, et al. Eur J Pharmacol, 2010, 629(1-3), 147-152.

Note: This product is for in vitro research use only