

AMG-458

Catalog No: #S6378


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	AMG-458
Brief Description	Inhibitors
Purification	/
Target Name	c-Met/HGFR inhibitor
Calculated MW	539.58
Formulation	C30H29N5O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: COC1=CC2=NC=CC(OC3=CN=C(NC(=O)C4=C(C)N(CC(C)(C)O)N(C4=O)C4=CC=CC=C4)C=C3)=C2C=C1
 [c:6,17,32,34,37,39,42,t:2,4,10,12,30] Pathways: Tyrosine Kinase/Adaptors Receptor: c-Met (D1228H); c-Met (H1094R); c-Met (Human); c-Met (Mouse); c-Met (V1092I) Boiling pt: Melting pt: Solubility: DMSO: 21 mg/mL (38.91 mM) Ethanol:

Note: This product is for in vitro research use only