

Amprenavir

Catalog No: #S6388


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amprenavir
Brief Description	Inhibitors
Purification	98.00%
Target Name	HIV Protease inhibitor
Calculated MW	505.63
Formulation	C25H35N3O6S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:CC(C)CN(C[C@@H](O)[C@H](CC1=CC=CC=C1)NC(=O)O[C@H]1CCOC1)S(=O)(=O)C1=CC=C(N)C=C1
 |r,c:12,14,35,t:10,30,32|Pathways:Microbiology&virology; Proteases/ProteasomeReceptor:HIV proteaseBoling pt:Melting pt:Solubility:DMSO:16 mg/mL
 (31.64 mM)Ethanol:100 mg/mL (197.77 mM)Water:

References

1. Helsley RN, et al. Mol Pharmacol, 2013, 83(6), 1190-1199.

Note: This product is for in vitro research use only