


Description

Product Name	AS-604850
Brief Description	Inhibitors
Purification	99.00%
Target Name	PI3K inhibitor
Calculated MW	285.22
Formulation	C11H5F2NO4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
SMILES: C1=CC2=C(C=C1/C=C\3C(=O)NC(=O)S3)OC(O2)F
Pathways: PI3K/Akt/mTOR signaling
Receptor: PI3K α ; PI3K β ; PI3K γ ; PI3K δ
Boiling pt: Melting pt: Solubility: DMSO: 57 mg/mL (199.84 mM) Ethanol: 5 mg/mL (17.53 mM) Water:

References

1. Camps M, et al, Nat Med, 2005, 11(9), 936-943.

Note: This product is for in vitro research use only