

Azaperone


Catalog No: #S6398

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Azaperone
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	327.4
Formulation	C19H22FN3O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:FC1=CC=C(C=C1)C(=O)CCCN1CCN(CC1)C1=CC=CC=N1 |c:3,5,22,24,t:1,20|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:65 mg/mL (198.53 mM)Ethanol:34 mg/mL (103.84 mM)Water:

Note: This product is for in vitro research use only