

Bepotastine Besilate

Catalog No: #S6405


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bepotastine Besilate
Brief Description	Inhibitors
Purification	99.83%
Target Name	Histamine Receptor antagonist
Calculated MW	547.06
Formulation	C ₂₁ H ₂₅ CIN ₂ O ₃ C ₆ H ₆ O ₃ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C1CN(CCC1O[C@@H](C2=CC=C(C=C2)Cl)C3=CC=CC=N3)CCCC(=O)O.C1=CC=C(C=C1)S(=O)(=O)OPathways: GPCR/
G Protein; Immunology/Inflammation; Neuroscience Receptor: H1 receptor Boiling pt: Melting pt: Solubility: DMSO: 109 mg/mL (199.24 mM) Ethanol:

References

1. Da Prada M, et al. J Pharmacol Exp Ther, 1989, 248(1), 400-414.

Note: This product is for in vitro research use only