

Beta-Lapachone


Catalog No: #S6407

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-Lapachone
Brief Description	Inhibitors
Purification	98.00%
Target Name	IDO inhibitor
Calculated MW	242.27
Formulation	C15H14O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism
SMILES: CC1(C)CCC2=C(O1)C1=CC=CC=C1C(=O)C2=O
Pathways: Metabolism
Receptor: IDO1
Boiling pt: Melting pt:
Solubility: DMSO: 33 mg/mL warmed (136.21 mM) Ethanol: 10 mg/mL (41.27 mM) Water:

References

1. Li CJ, et al. J Biol Chem. 1993, 268(30), 22463-22468.

Note: This product is for in vitro research use only