

Bexarotene

Catalog No: #S6410


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bexarotene
Brief Description	Inhibitors
Purification	99.80%
Target Name	Retinoid Receptor inhibitor
Calculated MW	348
Formulation	C24H28O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

systemSMILES:CC1=CC2=C(C=C1C(=C)C1=CC=C(C=C1)C(O)=O)C(C)(C)CCC2(C)CPathways:MetabolismReceptor:RXRBoling pt:Melting pt:Solubility:DMSO: 8 mg/mL (22.98 mM)Appearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhang C, et al. Clin Y Res, 2002, 8(5), 1234-1240.

Note: This product is for in vitro research use only