

Biapenem

Catalog No: #S6411


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Biapenem
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	350.39
Formulation	C15H18N4O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:[H][C@@]12[C@@H]([C@@H](C)O)C(=O)N1C(C([O-])=O)=C(SC1CN3C=NC=[N+]3C1)[C@@H]2C
|r,c:19,21,t:13|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:

References

1. Aldridge KE, et al. Antimicrob Agents Chemother, 1994, 38(4), 889-893.

Note: This product is for in vitro research use only