

Bindarit

Catalog No: #S6413


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bindarit
Brief Description	Inhibitors
Purification	98.71%
Target Name	CCR inhibitor
Calculated MW	324.37
Formulation	C ₁₉ H ₂₀ N ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism SMILES: CC(C)(OCC1=NN(CC2=CC=CC=C2)C2=C1C=CC=C2)C(O)=O

|c:11,13,16,19,21,t:5,9|Pathways: Immunology/Inflammation; Microbiology&virology Receptor: MCP-1/CCL2; MCP-2/CCL8; MCP-3/CCL7

Bolting pt: Melting pt: Solubility: DMSO: 65 mg/mL (200.38 mM) Ethanol: 65 mg/mL (200.38 mM) Water:

References

1. Zoja C, et al. *Kidney Int*, 1998, 53(3), 726-734.

Note: This product is for in vitro research use only