

## Bisacodyl

Catalog No: #S6414


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Bisacodyl
Brief Description	Inhibitors
Purification	97.13%
Target Name	Others
Calculated MW	361.39
Formulation	C22H19NO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:CC(=O)OC1=CC=C(C=C1)C(C1=CC=C(OC(C)=O)C=C1)C1=CC=CC=N1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 72 mg/mL (199.23 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Ikarashi N, et al. Am J Physiol Gastrointest Liver Physiol, 2011, 301(5), G887-895.

Note: This product is for in vitro research use only