

Caffeic Acid Phenethyl Ester


Catalog No: #S6429

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caffeic Acid Phenethyl Ester
Brief Description	Inhibitors
Purification	98.00%
Target Name	NF-κB inhibitor
Calculated MW	284.31
Formulation	C17H16O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: C1=CC=C(C=C1)CCOC(=O)/C=C/C2=CC(=C(C=C2)O)O
Pathways: NF-κB Receptor; NF-κB
Boiling pt: Melting pt: Solubility: DMSO: 57 mg/mL (200.48 mM) Ethanol: 57 mg/mL (200.48 mM) Water:

References

1. Natarajan K, et al. Proc Natl Acad Sci U S A. 1996, 93(17), 9090-9095.

Note: This product is for in vitro research use only