

Cefoperazone

Catalog No: #S6437


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cefoperazone
Brief Description	Inhibitors
Purification	99.30%
Target Name	Others
Calculated MW	645.67
Formulation	C25H27N9O8S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:othersSMILES:[H][C@]12SCC(CSC3=NN=NN3C)=C(N1C(=O)[C@H]2NC(=O)[C@H](NC(=O)N1CCN(CC)C(=O)C1=O)C1=CC=C(O)C=C1)C(O)=O |r,c:9,13,44,t:7,39,41|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:10 mM in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only