

Cetylpyridinium Chloride


Catalog No: #S6439

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cetylpyridinium Chloride
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	339.99
Formulation	C21H38N.Cl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:[Cl-].CCCCCCCCCCCCCCCC[N+]1=CC=CC=C1 |c:18,20,t:16|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:68 mg/mL (200.0 mM)Ethanol:68 mg/mL (200.0 mM)Water:68 mg/mL (200.0 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only