

Closantel Sodium

Catalog No: #S6452


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Closantel Sodium
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	685.06
Formulation	C22H13Cl2I2N2O2.Na
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:[Na+].CC1=C(NC(=O)C2=CC(I)=CC(I)=C2[O-])C=C(Cl)C(=C1)C(C#N)C1=CC=C(Cl)C=C1

|c:1,9,12,19,30,t:6,16,25,27|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:100 mg/mL (145.97 mM)Ethanol:100 mg/mL (145.97 mM)Water:

Note: This product is for in vitro research use only