

Cysteamine HCl


Catalog No: #S6460

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cysteamine HCl
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	113.61
Formulation	C2H7NS.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:Cl.NCCSPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:23 mg/mL (202.44 mM)Ethanol:23 mg/mL (202.44 mM)Water:23 mg/mL (202.44 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only