

MRP-L54 Polyclonal Antibody

Catalog No: #41177

Package Size: #41177-1 50ul #41177-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	MRP-L54 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	MRP-L54 Polyclonal Antibody detects endogenous levels of MRP-L54 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MRP-L54.
Target Name	MRP-L54
Other Names	MRPL54; 39S ribosomal protein L54; mitochondrial; L54mt; MRP-L54
Accession No.	Swiss-Prot: Q6P161NCBI Gene ID: 116541
Uniprot	Q6P161
GeneID	116541;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using MRP-L54 Polyclonal Antibody

Note: This product is for in vitro research use only