

Deltarasin

Catalog No: #S6462


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Deltarasin
Brief Description	Inhibitors
Purification	/
Target Name	PDE inhibitor
Calculated MW	603.75
Formulation	C40H37N5O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:C(OC1=CC=C(C=C1)C1=NC2=C(C=CC=C2)N1CC1=CC=CC=C1)[C@H](C1CCNCC1)N1C2=CC=CC=C2N=C1C1=CC=CC=C1
 |r,c:4,6,13,15,23,25,39,41,44,49,51,t:2,9,11,21,37,47|Pathways:MetabolismReceptor:PDEδBoiling pt:Melting pt:Solubility:DMSO:100 mg/mL (165.63 mM)Ethanol:100 mg/mL (165.63 mM)Water:

Note: This product is for in vitro research use only