

## Dexmedetomidine HCl

Catalog No: #S6466


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dexmedetomidine HCl
Brief Description	Inhibitors
Purification	99.86%
Target Name	Adrenergic Receptor agonist
Calculated MW	236.74
Formulation	C13H16N2.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:psychosisSMILES:Cl.C[C@H](C1=CN=CN1)C1=CC=CC(C)=C1C |r,c:4,10,13,t:2,8|Pathways:GPCR/G

ProteinReceptor:α2-adrenoceptorBoiling pt:Melting pt:Solubility:DMSO:48 mg/mL (202.75 mM)Ethanol:48 mg/mL (202.75 mM)Water:

## References

1.Xue F, Zhang W, Chu HC.Kidney Int. 2016 May;89(5):1164.

Note: This product is for in vitro research use only