

Diclofenac Potassium


Catalog No: #S6470

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Diclofenac Potassium
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	334.24
Formulation	C ₁₄ H ₁₀ Cl ₂ KNO ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: [K+].[O-]C(=O)CC1=CC=CC=C1NC1=C(Cl)C=CC=C1Cl | c:6,8,12,15,17,t:4 | Pathways: Others Receptor: Others Boling
 pt: Melting pt: Solubility: DMSO: 69 mg/mL (206.43 mM) Ethanol:

Note: This product is for in vitro research use only