

Difluprednate

Catalog No: #S6472


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Difluprednate
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	508.55
Formulation	C27H34F2O7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: [H][C@@]12CC[C@](OC(=O)CCC)(C(=O)COC(C)=O)[C@@]1(C)C[C@H](O)[C@@]1(F)[C@@]2([H])C[C@H](F)C2=CC(=O)C=C[C@]12C |r,c:36,t:32| Pathways: Others Receptor: Others Boling pt: Melting pt: Solubility: DMSO: 102 mg/mL (200.57 mM) Ethanol: 25 mg/mL (49.15 mM) Water:

Note: This product is for in vitro research use only