

Dioscin

Catalog No: #S6473


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dioscin
Brief Description	Inhibitors
Purification	99.85%
Target Name	Others
Calculated MW	869.04
Formulation	C45H72O16
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@@H]1CC[C@@]2([C@H]([C@H]3[C@@H](O2)C[C@@H]4[C@@]3(CC[C@H]5[C@H]4CC=C6[C@@]5(CC[C@@H](C6)O[C@H]7[C@@H]([C@H]([C@@H]([C@H](O7)CO)O[C@H]8[C@@H]([C@@H]([C@H]([C@@H](O8)C)O)O)O)O[C@H]9[C@@H]([C@@H]([C@H]([C@@H](O9)C)O)O)C)C)OC1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 22 mg/mL (25.31 mM); DMSO: 100 mg/mL (115.06 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang Z, et al. Biol Pharm Bull, 2001, 24(2), 159-162.

Note: This product is for in vitro research use only