

Divalproex Sodium


Catalog No: #S6474

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Divalproex Sodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Autophagy activator
Calculated MW	310.41
Formulation	C ₈ H ₁₆ O ₂ .C ₈ H ₁₅ O ₂ .Na
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: [Na+].CCCC(CCC)C(O)=O.CCCC(CCC)C([O-])=O
 Pathways: Autophagy
 Receptor: Autophagy
 Boling pt: Melting pt:
 Solubility: DMSO: 62 mg/mL (199.73 mM) Ethanol: 62 mg/mL (199.73 mM) Water: 62 mg/mL (199.73 mM)
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only