

DL-Carnitine HCl


Catalog No: #S6475

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DL-Carnitine HCl
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	197.66
Formulation	C7H15NO3.HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: [Cl-].C[N+](C)(C)CC(O)CC(O)=O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 40 mg/mL (202.36 mM)
Ethanol:

Note: This product is for in vitro research use only