

Domiphen Bromide

Catalog No: #S6477


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Domiphen Bromide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	414.46
Formulation	C ₂₂ H ₄₀ NOBr
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: [Br-].CCCCCCCCCCCC[N+](C)(C)CCOC1=CC=CC=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 83 mg/mL (200.26 mM) Ethanol: 83 mg/mL (200.26 mM) Water: 83 mg/mL (200.26 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Gao P, et.al. Eur J Pharmacol, 2012, 694(1-3), 45-52.

Note: This product is for in vitro research use only