

D-Pantothenic acid


Catalog No: #S6480

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D-Pantothenic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	219.24
Formulation	C9H17NO5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC(C)(CO)[C@H](C(=O)NCCC(=O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:44 mg/mL (200.69 mM)Ethanol:

References

1. Slyshenkov VS, et al. FEBS Lett, 2004, 569(1-3), 169-172.

Note: This product is for in vitro research use only