

Droxinostat


Catalog No: #S6481

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Droxinostat
Brief Description	Inhibitors
Purification	/
Target Name	HDAC inhibitor
Calculated MW	243.69
Formulation	C ₁₁ H ₁₄ ClNO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
SMILES: CC1=CC(Cl)=CC=C1OCCCC(=O)NO | c:4,6,t:1 | Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair;
NF- κ B Receptor: HDAC10; HDAC3; HDAC6; HDAC8; HDAC9
Boiling pt: Melting pt: Solubility: DMSO: 49 mg/mL (201.07 mM) Ethanol: 49 mg/mL (201.07 mM) Water:

Note: This product is for in vitro research use only