

Efaproxiral Sodium

Catalog No: #S6482


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Efaproxiral Sodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	363.38
Formulation	C ₂₀ H ₂₃ NO ₄ .Na
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:[Na+].CC1=CC(NC(=O)CC2=CC=C(OC(C)(C)C([O-])=O)C=C2)=CC(C)=C1

|c:19,21,24,t:1,8,10|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:73 mg/mL (200.89 mM)Ethanol:73 mg/mL (200.89 mM)Water:73 mg/mL (200.89 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Donnelly ET, et al. Exp Biol Med (Maywood), 2006, 231(3), 317-321.

Note: This product is for in vitro research use only