

Encorafenib (LGX818)

Catalog No: #S6487


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Encorafenib (LGX818)
Brief Description	Inhibitors
Purification	/
Target Name	Raf inhibitor
Calculated MW	540.01
Formulation	C22H27ClFN7O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: C[C@@H](CNC1=NC=CC(=N1)C2=CN(N=C2C3=CC(=CC(=C3F)NS(=O)(=O)C)Cl)C(C)C)NC(=O)OC Pathways: MAPK

Signaling Receptor: B-Raf (V600E) Boiling pt: Melting pt: Solubility: DMSO: 100 mg/mL (185.18 mM) Ethanol: 100 mg/mL (185.18 mM) Water:

Note: This product is for in vitro research use only