

Epinastine HCl


Catalog No: #S6488

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Epinastine HCl
Brief Description	Inhibitors
Purification	/
Target Name	5-HT Receptor inhibitor
Calculated MW	285.77
Formulation	C16H15N3.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:Cl[H].NC1=NCC2N1C1=CC=CC=C1CC1=CC=CC=C21 |c:10,12,18,t:2,8,16,20|Pathways:GPCR/G Protein;
NeuroscienceReceptor:HTBoling pt:Melting pt:Solubility:DMSO:57 mg/mL (199.46 mM)Ethanol:57 mg/mL (199.46 mM)Water:57 mg/mL (199.46 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only