

## Epinephrine HCl


Catalog No: #S6489

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Epinephrine HCl
Brief Description	Inhibitors
Purification	/
Target Name	Adrenergic Receptor inhibitor
Calculated MW	219.67
Formulation	C <sub>9</sub> H <sub>13</sub> NO <sub>3</sub> .HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system SMILES: Cl.CNC[C@H](O)C1=CC=C(O)C(O)=C1 |r,c:11,t:5,7|Pathways: Neuroscience; GPCR/G

Protein Receptor: Adrenergic Receptor Boiling pt: Melting pt: Solubility: DMSO: 44 mg/mL (200.3 mM) Ethanol: 7 mg/mL (31.86 mM) Water: 44 mg/mL (200.3 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only